




Discover and Play with the ZE5™ Cell Analyzer



Bio-Rad ZE5™ has been designed taking into account Core Managers indications , to overcome limitations of existing flow cytometers. It comes with:

- “ **Up to 30 optical parameters** for maximum flexibility and high-color multiplexing
- “ **Ability to Switch from different supports** (tubes, plates...) with ease
- “ **High Speed acquisition** for rare events analysis
- “ **Automatic functions** for daily use convenience
- “ **Powerful but simple user interface**

<p>Where:</p>			<p>OCEANOMED Campus de Luminy 163 Avenue de Luminy 13288 Marseille</p>
<p>When:</p>	<p style="text-align: center;">13-15 June 2018</p>		
<p>Program :</p>	<p>13 June 10h:30-13h:30- Introduction: (Andrea Tradori, Bio-Rad) Wetlab session (bring your samples!)</p> <p>14 June 10h:30-13h:30- Wetlab session (bring your samples!) Wetlab session (bring your samples!)</p> <p>15 June 9:30-13:00 Wetlab session (bring your samples!)</p>		

Please confirm your presence to:

Chrystele Romero : 06 72 26 85 43 chrystele_romero@bio-rad.com
Andrea Tradori: 06 45 63 90 49 andrea_tradori@bio-rad.com

**Please check that your fluorochrome combination is matching
ZE5's Lasers/Filters configuration!**